# Health o meter®

## P r o f e s s i o n a l

**ASSEMBLY & OPERATING INSTRUCTIONS  
DIRECTIVES DE MONTAGE ET D'UTILISATION  
INSTRUCCIONES DE MONTAJE Y USO**

**ELECTRONIC CHAIR SCALE  
MODEL 594KL  
BASCULA ELECTRONICA CON SILLA  
MODELO 594KL  
FAUTEUIL DE PESÉE ÉLECTRONIQUE  
MODÈLE 594KL**

Please register your product at [www.homscales.com](http://www.homscales.com)  
Registre su bascula en nuestro sitio web [www.homscales.com](http://www.homscales.com)  
Veuillez enregistrer votre produit sur [www.homscales.com](http://www.homscales.com)

**Thank you for your purchase of this product.** Please read this manual carefully and keep it handy for ready reference.

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## TABLE OF CONTENTS

CAUTION AND WARNING .....	2
SCALE SPECIFICATIONS .....	2
SET UP .....	3
OPERATION INSTRUCTIONS .....	3
PRESERVING BATTERY/PROGRAMMING AUTO-OFF TIME ..	3
MAINTENANCE .....	3
TROUBLESHOOTING .....	4
CALIBRATION PROCEDURE .....	4
LIMITED WARRANTY .....	5

## CAUTION AND WARNING

**To prevent injury and damage to your scale, please follow these instructions very carefully**

- Do not transport the scale while the chair is loaded.
- For accurate weighing the scale must be placed on a flat, stable surface.
- For accurate weighing, verify proper operation according to the procedure described in this manual before each use.
- Do not use in the presence of flammable materials.
- Operating at other voltages and frequencies than specified could damage the equipment.
- Do not leave patient unattended in the chair.
- Do not use chair scale in place of a wheelchair to transport patients.
- Do not bounce scale down stairs, as doing so may damage internal parts.
- Do not use the footrest as a step.

## SCALE SPECIFICATIONS

Health o meter Professional's Electronic Chair Scale, model 594KL, uses sophisticated microprocessor technology. This precision instrument is designed to provide accurate, reliable and repeatable weight measurements and features that make the weighing process simple, fast and convenient. The scale is set up to determine the weight of a stationary patient. The weight can be displayed in pounds or in kilograms. The unit can be operated using its 6 Volt rechargeable battery.

### SCALE SPECIFICATIONS

Capacity and resolution	600 lb x 0.2 lb / 270 kg x 0.2kg
Power Requirements	6 Volt Rechargeable battery, charger and battery included
Environmental	Operating temperatures: 5°C to 40°C (41°F to 104°F) Storage temperature: -1°C to +52°C (30°F to 125°F) Humidity: 85% rH
Chair	18 1/4" (w) x 15" (d) x 17 1/2" (h)
Scale Dimensions	24 3/8" (w) x 32 1/2" (d) x 36 3/4" (h) (Arm Width 20 1/2")

# SET UP

**Your scale is a precision instrument and should be handled accordingly.** The Health o Meter Professional Electronic Chair Scale is shipped completely assembled, so no additional assembly required.

Some protective plastic, cardboard and foam will need to be removed before using.

1. Remove metal bands attached to small skid.
2. Remove cardboard at the front and the back end of the chair.
3. Remove protective foam and plastic cover on the display.
4. Plug the scale's battery charger into the chair scale display and into a power source to charge the battery. The battery needs to be charged overnight to gain full power. When charging, be sure the On/Off switch located on the side of the display is in the On "I" position.

# OPERATING INSTRUCTIONS

1. Make sure the scale is on a hard, level surface before using.
2. Press the ON / OFF button to turn on the scale.
3. The indicator will start the initialization count from 000000 to 999999.
4. Wait until "0.0" and the arrows above ZERO and STABLE appear on the display.
5. The model 594KL model is designed to provide convenient and accurate weighing to the less mobile. The side-arms and footrest are movable for easy access and patient comfort.
6. Place the patient on the scale. The scale will display the patient weight. (NOTE: the weight will change until the patient is stable). Once the patient is stable an arrow will appear above STABLE on the display.
7. Press and hold the KG / LB button for 3 seconds to select the weigh mode, either kilograms or pounds.
8. Weight may remain displayed on the scale by pressing the HOLD / RELEASE button. To release, press the HOLD / RELEASE button again.

NO.	DESCRIPTION	FUNCTION
1.	HOLD / RELEASE	Holds the value of the weighed object on the display until the button is pressed again to
2.	BACK LIGHT	Turns display back light on and off
3.	KG / LB	Pressing and holding button for 3 seconds toggles between kilograms or pounds.
4.	TARE / ZERO	Zeros the scale prior to weighing or sets a tare value.
5.	ON / OFF	Turns scale ON and OFF

# PRESERVING BATTERY/PROGRAMMING AUTO-OFF

The 594KL is equipped with a power switch located on the side of the display.  When the switch is moved to the On "I" position it is supplying power to the scale, even when the display head is off. To charge the battery this switch must be in the On "I" position. To preserve power, keep the the switch in the Off "O" position.

The 594KL also has an Auto-Off feature that will automatically turn the scale off after a period of inactivity. This battery saving feature is programmable by the user. The 594KL can be programmed to automatically power off from 0 (disabled) to 21 minutes. The default shut off value set by the factory is 2 minutes. Once the auto-off time is changed it cannot be changed back to the 2 minute default setting.

1. Press ON / OFF and wait for the display to show "0.0".
2. At the same time press both the TARE/ZERO and KG/LB buttons. The display will show "t" along with a number representing the current shut-off setting.
3. Use the KG/LB button to change the shut off setting. The shut-off values are set in 3 minute increments. ("t 0" = disabled; "t 1" = 3 minutes; "t 2" = 6 minutes ..... "t 7" = 21 minutes.)
4. Press TARE/ZERO to save the desired turn off value. "CAL" will appear on the display.
5. Press BACK LIGHT to exit. Wait for "0.0" to appear and the scale is now ready for normal operation.

## GENERAL

The Maintenance and Troubleshooting sections provide instructions for maintenance, cleaning, troubleshooting and operator replaceable parts for the Electronic Chair Scale. Maintenance operations other than those described in these sections should be performed by qualified service personnel.

# MAINTENANCE

Before first-use and after periods of non-use, check the scale for proper operation and function. If the scale does not operate correctly, refer to qualified service personnel.

1. Check overall appearance of the total scale for any obvious damage, abuse, etc.
2. Inspect condition of battery charger.

## CLEANING

Proper care and cleaning is essential to ensure a long life of accurate and effective operation. Disconnect the scale from the battery charger.

1. Clean all external surfaces with a clean damp cloth or tissue. Mild soap and water solution may be used. Dry with a clean soft cloth.
2. Do not immerse the scale into cleaning or other liquid solution.
3. Do not use Isopropyl Alcohol or other solutions to clean the display surface.

## TROUBLESHOOTING

Refer to the following instructions to check and correct any failure before contacting service personnel.

SYMPTOM	POSSIBLE CAUSE	CORRECTIVE ACTION
Scale does not turn on	1. Dead battery 2. Faulty electrical outlet 3. Bad battery charger	1. Attach battery charger to recharge battery. 2. Use a different outlet for the battery charger. 3. Replace battery charger.
Questionable weight or the scale does not zero	1. External object interfering with the scale. 2. The display did not show "0.0" before weighing 3. Scale is not on a level floor 4. Scale is out of calibration	1. Remove interfering object from the scale. 2. Ask the patient to get off the scale and begin weighing process again. 3. Place the scale on a level floor and begin weighing process again. 4. Check weight with known weight value.
The display shows "- 0F- "	The load on the scale exceeds the maximum capacity (600 lbs / 270 kg)	Remove the excess weight and use the scale according to its limits.

## ERROR LIST

ERROR NO.	DEFINITION	CHECK / REPLACE
Err 1	A / D value is lower when CAL	Scale is below minimum calibration weight. Calibration weight must be more than 44 lbs. The recommended calibration weight is 600 lbs.
--OF--	The load on the scale exceeds the capacity (600 lb / 270kg).	Remove excess weight and use the scale according to its limits.
Or-Err	High saturation state (High A / D)	The load cell is not connected properly. Check cables and mechanical connection. If problem persists, contact Health o meter Professional Customer Service 800-815-6615.
Or-Err	Low saturation state (Low A / D)	The load cell is not connected properly. Check cables and mechanical connection. If problem persists, contact Health o meter Professional Customer Service 800-815-6615.
UF	Battery is low.	Recharge battery. Overnight charge is required to fully charge battery.

## CALIBRATION PROCEDURE

**This scale has been calibrated at the factory and does not require re-calibration prior to use.**

The following procedure should be performed only to correct any error in calibration. For maximum accuracy, calibration should be done using certified weights. The recommended calibration weight is 600 lbs, however calibration can be performed using any weight no less than 44 lbs and no more than 600 lbs.

1. Turn the scale on.
2. At the same time press the TARE/ZERO and KG/LB buttons. "t" and the number representing the current automatic shut off value will appear on the display.
3. Press the TARE/ZERO button. "CAL" and then '600.0' will appear on the display.
4. To change the calibration weight, press the HOLD/RELEASE button to move the flashing indicator to the left. Then press the KG/LB button to set the value to the desired calibration weight.
5. Press the TARE/ZERO button again. The display will show 'CAL 00'.
6. Press the TARE/ZERO button again. The display will show CAL-SP. Load 600 LB (or altered value) of calibration weight on to the scale and press the TARE/ZERO button. The display will show "600.0" (or altered value).
7. Remove the weight from the scale. The scale returns to the normal weighing mode.

# LIMITED WARRANTY

## What does the Warranty Cover?

Pelstar LLC scales are warranted from date of purchase against defects of materials or in workmanship for a period of one (2) years. If product fails to function properly, return the product, freight prepaid and properly packed to Pelstar. See "To Get Warranty Service" below for instructions. If manufacturer determines that a defect of material or in workmanship exists, customers' sole remedy will be repair or replacement of scale at no charge. Replacement will be made with a new or remanufactured product or component. If the product is no longer available, replacement may be made with a similar product of equal or greater value. All replaced parts are covered only for the original warranty period.

## Who is Covered?

The original purchaser of the product must have proof of purchase to receive warranty service. Please save your invoice or receipt. Pelstar dealers or retail stores selling Pelstar products do not have the right to alter, or modify or any way change the terms and conditions of this warranty.

## What is Excluded?

Your warranty does not cover normal wear of parts or damage resulting from any of the following: negligent use or misuse of the product, use on improper voltage or current, use contrary to the operating instruction, abuse including tampering, damage in transit, or unauthorized repair or alternations. Further, the warranty does not cover natural disasters, such as fire, floods, hurricanes and tornadoes. This warranty gives you specific legal rights, and you may also have other rights that vary from country to country, state to state, province to province or jurisdiction to jurisdiction.

To get warranty service make sure you keep your sales receipt or document showing proof of purchase. Call 1 (800) 638-3722 or 1 (708) 377-0600 to receive a return authorization number. Attach proof of purchase to your defective product along with your name, address, daytime telephone number and description of the problem. Carefully package the product and send with shipping and insurance prepaid to:

Pelstar LLC  
Attention R / A# \_\_\_\_\_  
Return Department  
9500 W. 55<sup>th</sup> St.  
McCook, IL 60525

If your scale is not covered by warranty, or has been damaged, an estimate of replacement costs will be provided to you for approval prior to replacing.

**PELSTAR, LLC**  
**9500 W. 55<sup>th</sup> St. • McCook, IL 60525 • USA**  
**1-800-638-3722 or 1-708-377-0600**

**PLEASE REGISTER YOUR SCALE FOR WARRANTY COVERAGE AT:**  
**[www.homscales.com](http://www.homscales.com)**

Health o meter® is a registered trademark of Sunbeam Products, Inc. used under license  
Health o meter® Professional products are manufactured, designed and owned by Pelstar, LLC.  
We reserve the right to improve, enhance or modify Health o meter® Professional product features or specifications without notice.

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**PELSTAR LLC**